


ScanPoint 3.0 image management technology provides viewing, printing and HIPAA-compliant online storage of patient exams, including images and voice annotations, generated from clinical use of AortaScan®, BladderScan®, and FloPoint® Elite instruments.

ScanPoint 3.0 lets you link your patients, exams, and electronic health record (EHR) systems like never before, through a new, improved, and easy-to-use screen interface.



Type	Model	Date/Time	Measurement
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:24:24 AM	168 ml
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:24:06 AM	171 ml
FloPoint® Elite	FloPoint® Elite S/N: 00001222	3/11/2010 10:23:16 AM	20.6 ml/s
AortaScan®	BV1 9600 S/N: 00002094	3/11/2010 10:15:18 AM	<3.0 cm
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:14:48 AM	0 ml
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:14:25 AM	0 ml
AortaScan®	BV1 9600 S/N: 00002094	3/11/2010 10:14:09 AM	<3.0 cm
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:13:46 AM	0 ml
AortaScan®	BV1 9600 S/N: 00002094	3/11/2010 10:01:49 AM	3.2 cm
BladderScan®	BV1 9600 S/N: 00002094	3/11/2010 10:01:05 AM	161 ml

First Previous 1 2 3 4 5 6 7 Next Last

Typical ScanPoint 3.0 exam screen

### ScanPoint 3.0 Benefits:


- Review, annotate, print, and archive exam data (images and voice annotations)
- Back up exam data automatically; no need to schedule or manage
- Securely interface with EHR systems to integrate exam data into patient health records
- Optimize instrument performance with online calibration and software upgrades
- Track calibration and warranty (service plan) status of each ScanPoint-enabled instrument



### ScanPoint 3.0 Supports:

- AortaScan AMI 9700 aortic measurement instruments
- BladderScan BVI 6000 series and BVI 9000 series bladder volume instruments
- FloPoint Elite uroflow instruments





Typical ScanPoint 3.0 BladderScan® Results screen

**ScanPoint 3.0 Features:**

- Supports HIPAA-compliant, secure access to exam records
- Provides a comprehensive view of individual patient exam results from multiple ScanPoint-enabled instruments, allowing physicians to quickly and easily track progress of treatment
- Generates PDFs of exam records in multiple formats for patient records or reimbursement
- Can be configured to notify biomedical technicians via e-mail of upcoming instrument calibration dates
- Fully localized user interface available in Chinese, English, French, German, Italian, Portuguese, and Spanish



Calibrate your BladderScan instrument online

**ScanPoint Service Plans**

ScanPoint service plans can be tailored to suit your needs. Your ScanPoint plan may include:

- ScanPoint Online Service – includes online archiving of patient exam records, EHR connectivity, Verathon® Customer Care support and online calibration and software upgrades.
- ScanPoint Online Service plus warranty coverage

A Premium Total Customer Care<sup>SM</sup> warranty is a cost-effective way to maintain your instrument. In addition to benefits of the ScanPoint online service, the premium warranty provides:

- a free Next Day loaner unit
- coverage of all related shipping charges within the USA
- coverage of repairs for accidental damage (a deductible may apply)
- instrument repair/replacement at no additional cost as determined by Verathon Technical Services

Please contact your local representative or Verathon Customer Care for pricing and further information.

**ScanPoint 3.0 Technical Specifications**

Compatible with Windows® XP operating system.



US Patent No. 6,569,097, 6,884,217. Other Patents Pending.  
ScanPoint, the ScanPoint symbol, AortaScan, the AortaScan symbol, BladderScan, the BladderScan symbol, FloPoint, Verathon, and the Verathon Torch are trademarks of Verathon Inc.  
© 2010 Verathon Inc. 0900-1083-03-86



**Corporate Headquarters**  
20001 North Creek Parkway / Bothell, WA 98011 USA  
425.867.1348 / 800.331.2313 / 425.883.2896 Fax  
[verathon.com](http://verathon.com)