

The GlideScope Direct intubation trainer is designed to facilitate instruction of classic direct laryngoscopy (DL).

As part of the GlideScope AVL family, the reusable Direct combines the characteristics of a standard Macintosh blade with AVL video technology. A digital video camera near the end of the blade and the digital AVL monitor allow instructors to watch and guide a direct laryngoscopic intubation as it happens.

The student sees an airway view as enabled by their DL technique, while the instructor observes the laryngeal view on the GlideScope AVL monitor.

The GlideScope Direct may be used to perform video-guided DL in a clinical setting, as well.

- ▶ Available in reusable and single-use configurations
- ▶ Digital camera and monitor for airway views in DVD clarity
- ▶ Real-time recording to help confirm tube placement and review DL intubation technique
- ▶ Advanced resolution output for OR theater viewing
- ▶ Anti-fogging mechanism to resist lens clouding and secretions
- ▶ An analog version of the Direct intubation trainer is also available for the GlideScope PGVL system



Now with single-use Stat option



## Designed to help instructors guide DL intubations.

### Specifications

Length (blade tip to handle):	119 mm
Height at camera:	15 mm
Width at camera:	23 mm
Max blade width in front of camera:	23 mm
Weight:	~600 grams

### Ordering P/N

Digital (for AVL):	0574-0103
Analog (for PGVL):	0574-0109
Direct Stat (box of 10):	0270-0712