



 GLIDESCOPE® Go™

 VERATHON®



# GO ANYWHERE.

PORTABILITY WITHOUT  
COMPROMISE.

[GlideScopeGo.com](http://GlideScopeGo.com)

# GO WITH CONFIDENCE.

Get the confidence of GlideScope® whenever and wherever you need it, with the portable GlideScope Go™ – the handheld system created for both EMS and mobile care teams.



**Intubate routine and difficult airways with the handheld system designed for rapid, first-pass success.**

## GO WITH CLARITY.

Get more of the light you want and less of light you don't. GlideScope Spectrum™ single-use video laryngoscopes deliver a high-resolution image powered by the latest advancements in lighting and camera technology, including Dynamic Light Control® and Ambient Light Reduction®, for enhanced airway views.

## GO WITH CONVENIENCE.

Compact, durable, easy to charge and ready on demand, GlideScope Go is ideal for use in small spaces, emergent procedures and whenever your environment demands mobility.

**Learn more at  
[GlideScopeGo.com](http://GlideScopeGo.com)**



## GO WITH COMPATIBILITY.

Access the portfolio\* of GlideScope® Spectrum™ single-use video laryngoscopes. Spectrum generates clear airway views and provides the widest range of single-use video laryngoscopes. Each is hot-swappable and compatible with your GlideScope cart-based or handheld system.



## GO WITH SMART FEATURES.

- Fully submersible video monitor and battery to simplify cleaning
- Easy to use and easy to train, with on/off button and intuitive display
- Rugged design ensures functionality even when dropped from 4 feet
- IP67 rating ensures full dust protection and reliable operation in wet conditions

- 3.5" landscape, color display with vertical tilt adjustment
- Scratch-resistant screen with anti-glare coating for reliable clarity in any environment

- Integrated battery charges to 80% in 2 hours and delivers a minimum of 100 minutes of continuous use on a full charge
- Operates on corded, external power should battery be discharged
- Remaining Charge indicator eliminates guessing
- Automatic recording with removable micro-USB drive
- Backed by a full warranty and customer support



\*See Operations and Maintenance Manual for full line of Spectrum video laryngoscopes that can be used with GlideScope Go.



## GO WITH THE SYSTEM YOU TRUST.

From the trusted leader in emergent and procedural airway management, GlideScope® Go™ serves as the go-to system for EMS and mobile care teams. An essential complement to GlideScope cart-based systems and compatible with our wide range of Spectrum™ single-use video laryngoscopes.



GlideScope® Core™



GlideScope® Spectrum™  
Single-Use Video Laryngoscopes



GlideScope® Go™  
Handheld System



Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA and Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.  
Willem Fenengastraat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201