⊜BladderScan :10[™]

CERTAINTY IN YOUR INVESTMENT

BladderScan i10 Value Adds Up As Costs Go Down

BladderScan i10TM provides more than trusted accuracy. It also offers value-added features, service and warranty options to reduce your total cost of ownership.



INTRODUCTORY PRICE vs. TOTAL COST OF OWNERSHIP


VS

BladderScan i10

BladderScan i10 reduces cost of ownership over time with built-in reliability, durable industrial design and comprehensive warranty solutions.

PURCHASE PRICE BENEFITS INCLUDE:

- Comprehensive 5-year warranty
- No routine calibration needed
- CaliScan[™] probe self-diagnostics
- On-site and digital training resources
- Complimentary software updates
- Simple charting tool
- Staging and deployment tools




Competitor

If a low purchase price seems too good to be true, it probably is. Unexpected maintenance and frequent downtime can make your total cost of ownership a steep price to pay.

PURCHASE PRICE MAY NOT INCLUDE:

- Probe drops and replacements
- Calibration
- Training
- · Service and repair
- Preventive maintenance
- Software updates and licenses
- Loaners







BladderScan i10[™] and Verathon® – Working For You

It takes more than reliable product design to ensure around-the-clock uptime. With Verathon, full line technical support and customer service teams augment your on-site support. Complimentary on-site training, software updates and fleet management services optimize usage. You can keep out-of-pocket costs to a minimum over the life of your device, thanks to a comprehensive warranty covering all damage and failures associated with the probe and console. You can have peace of mind knowing your investment is working every day to improve patient care.

One Call Does It All

Whatever you need, our experienced local customer service and full-line technical support teams are ready to help. Plus, we provide same-day loaner shipping when your order is processed by 2:00 PM PST.



THE VALUE INCLUDED IN THE PRICE



TRUSTED ACCURACY FOR ANY USER

Powered by ImageSense[™], Verathon's groundbreaking artificial intelligence technology, BladderScan i10 empowers users to consistently scan with an accuracy of ± 7.5% on volumes greater than 100 mL (or ± 7.5 mL on volumes less than 100 mL)* across all adult and pediatric patient types and anatomies. The large touchscreen interface uses intuitive graphics, automatic bladder location detection and BladderTraq[™] contour tracking to guide every step.



ACCELERATED WORKFLOW

Optimized for roaming through hectic hospital corridors, BladderScan i10 pairs with a quiet glide workstation that features a 360-degree handle for easy positioning. Single-handed probe docking keeps users focused on the patient – not on managing the equipment. ImageSense, combined with VMODE® probe technology, automatically measures the bladder volume within seconds. With a simple charting tool, you can easily transfer exam results to any electronic medical record system.



TOUGH AND READY DEPENDABILITY

BladderScan i10 is designed to meet the new realities of infection control, using custom-blended medical grade materials to withstand fast-kill disinfectants. The rugged, fully-sealed IPX4 probe and dripresistant IPX2 console ensure your precision electronics are protected against accidental spills. No calibration is needed, even if you swap probes. CaliScan™ probe diagnostics confirm probe function in the event of a drop. BladderScan i10 also has both onboard charging and a user-swappable battery, so the system is always ready whenever and wherever you need it.

Lower your cost of ownership while improving patient care.

Contact your BladderScan® representative today or Verathon Customer Care at 800-331-2313.

*Assumes BladderScan i10 is being used in accordance with stated instructions, scanning a tissue-equivalent phantom.





Verathon Inc. 20001 North Creek Parkway Bothell, WA 98011 USA

Tel: +1 800 331 2313 (USA & Canada only)

Tel: +1 425 867 1348 Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V. Willem Fenengastraat 13 1096 BL Amsterdam The Netherlands

Tel: +31 (0) 20 210 30 91 Fax: +31 (0) 20 210 30 92 Verathon Medical (Australia) Pty Limited Unit 9, 39 Herbert Street St Leonards NSW 2065 Australia

Tel: 1800 613 603 (Within Australia) Fax: 1800 657 970

Tel: +61 2 9431 2000 (International) Fax: +61 2 9475 1201