

# Designed for Portability

The BladderScan® BVI 6100 bladder volume instrument, with VMODE® technology, is a portable 3D ultrasound device that quickly, accurately, and noninvasively measures urinary bladder volume.

## **Quick and Accurate**

measurements of urinary bladder volume and post-void residual (PVR)

# **Precision Aiming**

via lightweight, portable probe

# **Designed around you**

easy to use, distinct design

# **Handheld Portable Design**

The handheld BVI 6100 offers unprecedented portability and versatility in multiple clinical settings.





# **VMODE Technology**

In less than 5 seconds, VMODE® technology automatically captures 12 B-Mode slices of the bladder, and displays the calculated volume results. VMODE automated 3D quantification of bladder volume helps clinicians assess bladder function quickly and confidently.

# **Aiming Assistance**

Directional arrows on the probe screen assist the user in ensuring the most accurate measurement is taken.



### BladderScan.com



Verathon Inc. 20001 North Creek Parkway Bothell, WA 98011

Tel: +1 800 331 2313 (USA & Canada only)

Tel: +1 425 867 1348 Fax: +1 425 883 2896

#### EC REP

Verathon Medical (Europe) B.V. Willem Fenengastraat 13 1096 BL Amsterdam The Netherlands

Tel: +31 (0) 20 210 30 91 Fax: +31 (0) 20 210 30 92 Verathon Medical (Australia) Pty Limited Unit 9, 39 Herbert Street St Leonards NSW 2065 Australia

Tel: 1800 613 603 (Within Australia) Fax: 1800 657 970

Tel: +61 2 9431 2000 (International) Fax: +61 2 9475 1201