Trusted Accuracy Meets Tested Durability





BladderScan i10— The Next Generation Bladder Scanner

BladderScan i10[™], powered by ImageSense[®] technology, delivers trusted accuracy for every user, accelerated workflow, tough and ready dependability, and certainty in your investment.

A partnership with Verathon® goes beyond the device itself. We deliver a comprehensive range of services tailored to your needs, giving you the power to maximize BladderScan® usage, minimize your cost of ownership, and improve patient care.




Trusted Accuracy for Every User

Powered by ImageSense®, Verathon's groundbreaking artificial intelligence technology, BladderScan i10™ harnesses real-world data to consistently detect and define bladder location, size and shape across all pediatric and adult patient types and anatomies.

± 7.5% on volumes greater than 100 mL¹

± 7.5 mL on volumes less than 100 mL



CAUTI is largely preventable. Effective use of bladder scanner technology reduces unnecessary catheterization by nearly 50%², lowering the risk of infection.

Successful CAUTI reduction initiatives depend on broad use and adoption by patient care staff of reliable and dependable bladder scanner technology.



^{1.} Accuracy specifications assume BladderScan i10 is being used in accordance with stated instructions, scanning a tissue-equivalent phantom.

^{2.} Clinical Care Improvement Strategies, Preventing Catheter Associated Urinary Tract Infections, 2013 The Joint Commission.

Accelerated Workflow

ImageSense®, combined with VMODE® probe technology, automatically measures the bladder volume within seconds.

- Large 10.1" (25.6 cm) optically bonded touchscreen reduces glare
- Interface uses intuitive graphics, automatic bladder location detection and BladderTraq[™] contour tracking to guide the user through every step
- Simple charting tool easily transfers exam results to any electronic healthcare record system

Optimized for roaming through hectic hospital corridors, BladderScan i10 pairs with a quiet glide workstation.


- Single-handed probe docking keeps users focused on the patient – not on managing the equipment
- 360-degree handle enables easy positioning
- With both onboard charging and user-swappable battery, the system is always ready whenever and wherever you need it

Simple EHR Charting

Scan the touchscreen directly, or scan the printout.







Tough and Ready Dependability

BladderScan i10[™] is designed and built for continuous operation in taxing hospital environments.

- Durable construction consisting of medical grade materials, metal gears and industrial-strength locking cables ensure reliable operation and high service availability
- No routine calibration is needed, even if you swap probes
- CaliScan[™] probe diagnostics confirm probe function in the event of a drop
- BladderScan i10 passes a probe drop test that meets or exceeds a 4' (1.22m) drop to a steel surface using the MIL-STD 810H standard

BladderScan i10 is designed to meet the new realities of infection control.

- Custom-blended medical grade materials withstand fast-kill disinfectants
- Rugged, fully-sealed IPX4 probe and drip resistant IPX2 console ensure your precision electronics are protected

Clean with Ease

Chemically resistant plastics used at each critical touchpoint.





Product Highlights



Accuracy for Every User



Accelerated Workflow



Tough and Ready Dependability



Flexible Power Management



Click or scan to view 360°

Certainty in Your Investment

You can keep out-of-pocket costs to a minimum over the life of your device, thanks to a comprehensive factory warranty and exceptional support.

Complimentary on-site training, software updates and fleet management services optimize usage. Technical support and customer service teams are on call to augment your on-site support.

Verathon® enhances the value of BladderScan i10™ with exceptional support, so you can have peace of mind knowing your investment is working every day to improve patient care.



Comprehensive 5-Year Warranty

- 5-Year Factory Warranty covers all damage and failures associated with the probe and console
- Same-day loaner shipping when your order is processed by 2:00 PM PST
- "One call does it all". If you need assistance, our experienced US-based customer service and full-line technical support teams are ready to help.
- Typical service turnaround within 5-days from receipt of instrument

Optimized Uptime

Your customer experience is our top priority. Here are a few examples of how we're making your life easier:

- · No routine calibration needed
- Onboard diagnostics to support your preventative maintenance program
- Complimentary software upgrades allow you to update your instrument over time
- Quick transfer of configuration settings:
 Clone preferred device settings across the fleet
- Device configuration lockdown functionality ensures compliance with individual hospital regulations and best practices
- Service Partner Programs



"Verathon's on-site support has reinvigorated our bladder care program and helped us improve our CAUTI prevention practices."

Jonathan Kuehn, BSN, RN-BC Aspirus Wausau Hospital

Training and Education

It takes more than a great device to improve patient care. You also need the know-how that enables you to fully integrate the device into your practice. Verathon® provides the education and training to enable clinical adoption of your BladderScan® purchase:

- · On-site training
- Comprehensive training videos to educate on the product and bladder scanning skills
- In-service resources: posters, sign-in sheets, staff competency certification, and comprehensive overviews and quizzes

Let our team help you develop a plan that maximizes your BladderScan usage and streamline your operations so you can get the most out of your system.



Accessories





Barcode Scanner 0810-0291



Printer **0800-0640**



Smart Battery Charger **0400-0157**



Calibration Tank 0620-1083



Thermal Printer Paper (x 5) 0800-0319



Rechargeable Battery 0400-0155



Tissue Phantom 0620-0068



Ultrasound Gel 0800-0005



Ultrasound Gel 0800-0643 20g single use packette, Box of 100

Learn more about the full BladderScan® suite of products, clinical resources and In-service instructional videos at bladderscan.com



verathon.com



20001 North Creek Parkway Bothell, WA 98011 USA

Tel: +1 800 331 2313 (USA & Canada only)

Tel: +1 425 867 1348 Fax: +1 425 883 2896



Verathon Medical (Europe) B.V. Willem Fenengastraat 13 1096 BL Amsterdam The Netherlands

Tel: +31 (0) 20 210 30 91 Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited Unit 9, 39 Herbert Street St Leonards NSW 2065 Australia

Tel: 1800 613 603 (Within Australia)

Fax: 1800 657 970

Tel: +61 2 9431 2000 (International)

Fax: +61 2 9475 1201