




Designed to Work
with GlideScope
Video Laryngoscopes

GlideRite Stylets—Enabling Quick Intubations



Designed to Work Together

GlideRite® Stylets are designed to complement angulated GlideScope® video laryngoscopes and enable quick placement of ET Tubes.

- Rigid design enhances Endotracheal Tube control for fast intubation
- Preformed to match the shape of GlideScope angulated blades
- Reusable and single-use options
- Options available to cover a broad range of patients and procedures
- Available in packs of 10

Single-use Stylets

Packaged for immediate use.



GlideRite® Single-use Stylet – Small

Preformed, rigid biocompatible plastic

- Designed for use in ET Tubes 3.0 mm to 4.0 mm
- Adjustable stopper



GlideRite Single-use Stylet – Medium

Preformed, rigid stainless steel

- Designed for use in ET Tubes 4.5 to 5.5 mm
- Easy to use handle design for simple ET Tube securement and stylet removal



GlideRite Single-use Stylet – Large

Preformed, rigid stainless steel

- Designed for use in ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal

② Indicates single-use

Reusable Stylets

Compatible with high-level disinfectants and autoclave.



GlideRite Rigid Stylet

Preformed, rigid stainless steel

- Designed for use in ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal



GlideRite DLT Stylet

Preformed, rigid stainless steel

- Designed for use in double lumen ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal
- Unique handle design with integrated tracheal lumen positioning pin

Schedule a demo with your Verathon representative or visit [verathon.com/stylets](https://www.verathon.com/stylets)

verathon

[verathon.com](https://www.verathon.com)

Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA

Tel: +1 800 331 2313
(USA & Canada only)

Tel: +1 425 867 1348
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.
Willem Fenengastrat 13
1096 BL Amsterdam
The Netherlands

Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia

Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970

Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201