



# Portability without Compromise



# Clarity. Convenience. Ease.

GlideScope® Go™ 2 is Verathon's second generation handheld video laryngoscope system. It delivers the confidence, clarity and ease-of-use of GlideScope®, your trusted brand of video laryngoscopes designed for first pass success<sup>1,2</sup> in a wide range of patients.



## Go with Clarity

GlideScope Go 2 with Spectrum™ QC delivers a high-resolution image and enables clear airway views.

- Clear visualization on a scratch-resistant 3.5" Corning® Gorilla® Glass<sup>3</sup> screen
- Enhanced image quality with a wider field of view to see more



## Go with Convenience

Get the confidence of GlideScope wherever and whenever you need it. Compact and durable, GlideScope Go 2 moves with you.

- Accessible and ready on demand for routine and difficult airways
- Spectrum QC blades are interchangeable with GlideScope Core and Go 2 for standardization throughout your organization



## Go with Ease

GlideScope Go 2 includes intuitive, easy-to-use features along with the reliable GlideScope technology that you know and trust.

- Easily click together Go 2 with Spectrum QC or Video Baton QC Large with magnetic QuickConnect™ technology
- Customize your device settings using the intuitive touchscreen design



Learn More about  
Spectrum QC eco

## Innovation with the Environment in Mind

Looking for a sustainable option? Explore Spectrum QC eco™, made with an 80% bio-based plastic shell<sup>4</sup>.

1. Silverberg, Michael J., et al. "Comparison of Video Laryngoscopy Versus Direct Laryngoscopy During Urgent Endotracheal Intubation: A Randomized Controlled Trial." *Critical Care Medicine* 43.3 (2014): 636-641.  
2. Ruetzler K, Bustamante S, Schmidt MT, et al. Video Laryngoscopy vs Direct Laryngoscopy for Endotracheal Intubation in the Operating Room: A Cluster Randomized Clinical Trial. *JAMA*. 2024;331(15):1279–1286. doi:10.1001/jama.2024.0762  
3. Corning and Gorilla are registered trademarks of Corning Incorporated.  
4. The raw material of the shells is linked to 80% bio-circular feedstock, which is allocated via the mass balance approach based on the amount of bio-circular material sourced in its production.

# Your Go 2 Handheld Video Laryngoscope System

Compact and durable, GlideScope® Go™ 2 is ready when your environment demands mobility.

## Durable 3.5" Touchscreen

Go 2 is ready for your environment with Corning® Gorilla® Glass and IP67 rating

Color display on a landscape screen

## Charging Cradle

Simplify charging with optional charging cradle that can be mounted vertically

## Battery Indicator

Remaining charge is always visible

## Rechargeable Battery

Up to 100 minutes of continual use with the integrated battery

## Adjustable Hinge

Adjust viewing angle with 170° vertical tilt



## Magnetic QuickConnect

Enhances efficiency and connectivity

## Touchscreen Navigation

Enable settings such as auto-record, auto-shutdown, and Dynamic Light Control®

## Dynamic Light Control

Dynamic Light Control optimizes image quality and brightness at the vocal cords

## Onboard Storage

Capture snapshots and videos with storage for up to 40 hours of video

# Compatible with a Broad Portfolio of Video Laryngoscopes

Clear, convenient, and easy to use, GlideScope® Go™ 2 is compatible with Spectrum™ QC and Video Baton QC. It complements your GlideScope® Core™ advanced integrated system to deliver GlideScope on the Go.



## Spectrum QC

Compatible with Spectrum QC and Spectrum QC eco™.



Hyperangle S1



Hyperangle S2



Hyperangle S2.5



Hyperangle S3



Hyperangle S4



MAC S1



MAC S2



MAC S3



MAC S4



Miller S00



Miller S0



Miller S1



Miller S2



## Video Baton QC

Compatible with Video Baton QC and single-use GVL® Stats.



Video Baton QC Large



GVL 3 Stat



GVL 4 Stat

Spectrum QC and GVL Stat products are sold in boxes of 10.

Schedule a demo with your Verathon representative or visit [verathon.com/glidescopego2](http://verathon.com/glidescopego2)

**verathon**

verathon.com

Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA & Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.  
Willem Fenengstraat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201