


GLIDERITE STYLETS



ENABLING QUICK INTUBATIONS



DESIGNED TO WORK WITH GLIDESCOPE VIDEO LARYNGOSCOPES

The angle of GlideRite® Stylets complements the unique angle of GlideScope® video laryngoscopes to help facilitate quick placement of ET Tubes.

- Rigid design enhances Endotracheal Tube control for fast intubation
- Preformed to match GlideScope blade angle
- Reusable and single-use options available
- Available in packs of 10

TO LEARN MORE, CONTACT YOUR
LOCAL VERATHON REPRESENTATIVE
OR CALL 800-331-2313

SINGLE-USE STYLETS

Packaged sterile for immediate use



GLIDERITE® SINGLE-USE STYLET – SMALL

Preformed, rigid plastic

- Designed for use in ET Tubes 3.0 mm to 4.0 mm
- Adjustable stopper
- Bio-compatible plastic



GLIDERITE® SINGLE-USE STYLET – LARGE

Preformed, rigid stainless steel

- Designed for use in ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal

REUSABLE STYLETS


Compatible with high-level disinfectants and autoclave



GLIDERITE® RIGID STYLET

Preformed, rigid stainless steel

- Designed for use in ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal



GLIDERITE® DLT STYLET

Preformed, rigid stainless steel

- Designed for use in double lumen ET Tubes 6.0 mm and larger
- Easy to use handle design for simple ET Tube securement and stylet removal
- Unique handle design with integrated tracheal lumen positioning pin

② indicates single-use

glidescope.com



Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA
Tel: +1 800 331 2313
(USA and Canada only)
Tel: +1 425 867 1348
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.
Willem Fenengastrat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited
Unit 9, 39 Herbert Street
St Leonards NSW 2065
Australia
Tel: 1800 613 603 (Within Australia)
Fax: 1800 657 970
Tel: +61 2 9431 2000 (International)
Fax: +61 2 9475 1201